



Medical Electric Power Tools Model: BLJZ-1 Series 1000C Premium & Best Sale -Newest Ergonomic design, offering comfortable experience.

- -Environment friendly nickel-metal hydride battery.
- -Autoclavable up to 135°C at 1.6 MPa
- -High precise gear box, more stable performance
- -Imported motor, low noise, longer service life
- -Imported chuck

-New designed machine body to distinguish your brand



Introduction of products

1.1 Product formation

Medical electric saw drill is consisting of handpiece, battery and charger. The charger is not belonging to medical device. Warranty: 12 months for Handpiece; 3 months for battery.

1.2 Intended usage and characters

Electric saw drill is our newly developed product which is used for orthopedic medical institutions to drill or saw bone surgical operation. Characteristics of product are stepless speed regulation and at the level of waterproofing is IPX8.









Handpiece of a Drill Handpiece of a Saw Ba

Battery







Medical Electric Power Tools Model: BLJZ-1 Series 1000C Premium & Best Sale

Introduction of products

2 Technical parameters

2.1 Product list and basic parameter (please refer to next page)2.2 Charge

Input: AC 100~240V, 50/60Hz, 30VA Output: DC14.4V, 700mA;

2.3 Battery

2.3.1 Voltage: DC14.4V, Ni-Mh1800mAh2.3.2 Fully charged time: about 3 hours;2.3.3 Service time for full-load condition: ≥2minutes;

2.4 Environmental conditions Normal operation conditions

a)Ambient temperature range: $5^{\circ}C \sim 40^{\circ}C$;

b)Relative humidity range: ≤70%;

c)Atmospheric pressure range: 860hPa~1060hPa;

Transportation and storage condition

a)Ambient temperature range: $-10^{\circ}C \sim 40^{\circ}C$;

b)Relative humidity range: $\leq 90\%$;

c)Atmospheric pressure range: 500hPa~1060hPa;







Medical Electric Power Tools Model: BLJZ-1 Series 1000C Premium & Best Sale





Series 10	000 Series 10	000C Series 200	00 Series 5000	Series 8000	Series 8000
Part No.	BL1101C	BL1102C	BL1201C	BL1103AC	BL1204C
Picture					
Description	Cannulating Drill For Trauma Surgery	Oscillating Saw For Joint Surgery	Dual-function Cannulating Drill With quick-Change chuck For Trauma Surgery	Bone Drill For Trauma surgery	Dual-function Acetabulum Reaming Drill With quick-Change chuck For Joint Surgery
Parameter	No-load Speed(r.p.m) :700 Torque(Nm):≥5.0 Operation(V) voltage:14.4 Sterilization: 135°C Noise(db):≤75 Radial run-out:≤0.1mm	Oscillations(osc/min):≥12000 Swimgs(°):≥5 Operation(V) voltage:14.4 Sterrilization:135°C Noise(db):≤90	No-load Speed(r.p.m) :700 Torque(Nm):≥5.0 Operation(V) voltage:14.4 Sterilization: 135°C Noise(db):≤75 Radial run-out:≤0.2mm	No-load Speed(r.p.m) :1100 Torque(Nm):≥4.0 Operation(V) voltage:14.4 Sterilization: 135°C Noise(db):≤75 Radial run-out:≤0.1mm	No-load Speed(r.p.m) :450 Torque(Nm):≥10.0 Operation(V) voltage:14.4 Sterilization: 135°C Noise(db):≤75 Radial run-out:≤0.2mm

* Boly offers cannulated drill BL1101AC/BL1201AC with a speed of 1100rpm for customization, welcome to consult Boly for more details.



Series 1000

Series 1000C

Series 4000


Series 5000

Series 8000

Series 2000



BEST CONBINATION FOR JOINTS/HIP SURGERIES



www.bolymedical.com

Aluminum Carrying Case

Configurations:

Handpiece of BL1204C x1 Handpiece of BL1102C x1 Battery x2 Charger x1 Aseptic transfer kit x2 Drill Chuck x1 Key to chuck x1 Saw blades x3

Package:

Neutral Outer Carton, 1-case Contained (425*140*370mm)

Sterilization Case

Configurations:

(455x395x90mm)

Handpiece of BL1204C x1 Handpiece of BL1102C x1 Battery x2 Charger x1 Aseptic transfer kit x2 Drill Chuck x1 Key to chuck x1 Saw blades x3 Package: Neutral Outer Carton, 1-case Contained



Series 10	Series 10000 Accessories	Series 2000	Series 5000	Series 8000	Series 8000
Part No.	BL0001	BL0002	BL0003	BL0004A	BL0006A
Picture		Con Con			H a
Description	Nickel-metal hydride battery	Charger	Aseptic Transfer Kit	Кеу	Saw Blades
Quantity	2 pcs	1 pcs	1 pcs	1 pcs	3 pcs

